



# USER INSTRUCTIONS



## INTENDED USE

Masks can protect people from coming in contact with respiratory droplets and airborne particles, especially when exposed to coughing, sneezing or when near others who are speaking. The masks are effective even if an individual is infected with a virus or fatal pathogen but not showing symptoms, provided the mask is used correctly and consistently and is being properly worn and fitted for optimum safety and comfort.

## APPLICATION

### Protection for Oneself

Acts as a shield to reduce the wearer's risk from potential harm in connection with inhalation of airborne particles.

### Donning Mask




1 Before donning mask, wash your hands, use hand sanitizer, and/or put on a clean pair of gloves.



2 With the inside facing you and the nose wire on the top, place the loops around the ears.



3 Pull top and bottom of mask to expand the folds. The top of the mask should cover your nose, and the bottom should cover your chin.



4 Pinch or press the nose wire to it fit the shape of your nose bridge.

### Protection for Others

Acts as a containment barrier to reduce risk of contact with respiratory droplets and airborne particles.

1 Before doffing mask, wash your hands, use hand sanitizer and/or put on a clean pair of gloves.

2 Carefully remove mask from the face, touching only the loops.

3 Dispose of the mask in a waste container.

4 Wash your hands or use a hand sanitizer immediately after removing your mask



### Doffing Mask

**CAUTION:** The outside surface of the mask should never come in contact with your face; avoid pulling the mask up on your forehead and do not let it hang on one ear or around your neck or on your arm



Do not touch front of mask when wearing mask



Do not wear on nose



Do not wear on your chin